


Invacare® IVC™ Tracer® SX5 Recliner Wheelchair

Seat-to-Floor 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	Recliner Axle Mounting Bracket
18 1/3"	Middle	22"	Recliner Axle Mounting Bracket
19 1/2"	Bottom	24"	Recliner Axle Mounting Bracket

Standard on SX5 Recliner

- Rear anti-tippers
- Patient only operated push-to-lock wheel lock

Weight Limitation 250 lbs., 300 lbs. 20" & 22" width
 Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.



Comments

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation

USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2018 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.